

contacta 

# Single-bay Infrared Receiver/ Hearing Loop Listener Charger

IR-RX2-DC1

Installation &  
User Guide

**General information:**

Thank you for your purchase of the single-bay charger.

Please read through the operating instructions carefully to be able to set up the unit correctly and to prevent damage to the rechargeable batteries within the receiver.

**Included in the standard package:**

Please check if all of the following components are included:

- Single bay charger
- Power supply unit
- Manual

If any parts are missing, please immediately contact your dealer or the manufacturer directly.

**Function:**

This charger is able to recharge the batteries of a single bay infrared receiver/hearing loop listener charger.

**Setup**

Place the single bay charger in a dry and ventilated position. First connect the power supply unit with the charger, then plug it into a electrical mains outlet.

**Charging the receivers:**

When a receiver is placed into the charging bay, the control lamp lights up in red. After the receiver is fully charged, the control lamp will display a green light. It may take up to several hours until a receiver is fully charged.

The control lamp switches off when the receiver is removed from the bay.

**Warranty**

Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer or the manufacturer directly.

This warranty covers the repair of the product and returning it to you free of charge. The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people not authorised to do so (destruction of the seal on the unit).

Repairs will only be carried out under warranty if the completed warranty card is returned accompanied by a copy of the dealer's in voice/till receipt.

**Always specify the product number in any event.**



**Disposal** of used electric and electronic units (applicable in the countries of the European Union and other European countries with a separate collection system). The symbol on the product or the packaging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the recycling of electric and electronic units. You protect the environment and health of your fellow individuals by the correct disposal of this product. Environment and health are endangered by a faulty disposal. Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.

**Specifications:**

Power supply:	100-240 Volt AC 50/60 Hz
Number of charging bays:	
Control lamps:	Dual colour LED (red / green)
Dimensions:	L 85 mm x W 135 mm x H 32 mm
Weight:	approx. 190 g


This device complies with the following EU directives:



- Low Voltage Directive 2014/35/EC
- EMC Directive 2014/30/EC
- RoHS Directive 2011/65/EC
- WEEE Directive 2012/19/EC

Compliance with the directives listed above is confirmed by the CE seal on the device.

Technical specifications subject to change without prior notification.



Local dealer:

UK & ROW  
+44 (0) 1732 223900  
sales@contacta.co.uk

US & Canada  
+1 616 392 3400  
info@contactainc.com